



www.gmi-inc.com  
sales@gmi-inc.com  
(763) 712-8717



# Labculture® G4

## Class II Type B2 Biological Safety Cabinets

*The Most Advanced, Energy-efficient, Safe,  
and Ergonomic Biosafety Cabinet in the World*




# LABCULTURE® G4 (LB2 G4) CLASS II TYPE B2 CABINET,



## USB Port and Zero Volt Relay Contact

- USB Port to send operational information to Building Management System (BMS)
- Zero Volt Relay Contact to turn ON/OFF exhaust blower and signal the building alarm



## Airflow Sensor

- Monitors real-time airflow for safety
- Alert the user if airflow is insufficient



## Centurion 7" Capacitive Touchscreen Controller

- Displays all safety information on one large screen
- Shows cabinet parameters with intuitive 3D illustration
- Easy to use menu, similar to Smart Phone Apps
- Large buttons, easy to operate when wearing gloves
- Self-guidance to users to deal with specific situations
- Centered and angled down for easy reach and viewing
- Optional: 21 CFR Part 11 Compliance



## Single Piece Wall

- Easy to reach service fixtures and electrical outlets on sidewalls
- Large radius corners for easy cleaning



## User-friendly Work Tray

- Largest useable area in the market
- Recessed to contain spillage
- Sloped perimeter for easy cleaning
- Large, easy to clean tray handle



## Raised Arm Rest

- Prevent grille blocking
- Comfortable working posture
- Durable stainless steel construction



## Ergonomic Work Zone

- 10° angle to optimize user comfort, reduce glare, and maximize reach into the work area
- Brightly illuminated with >1200 lux (111 ft. cd)
- Industry-leading dimmable LED for optimum work comfort
- Airtight seal port for cable/tube exit protected by a negative pressure side wall





Esco Labculture® G4 Class II Type B2 Biosafety Cabinet  
Available in 3 feet, 4 feet, 5 feet, and 6 feet models.

# FEATURING ADVANCED TOUCHSCREEN CONTROLLER

## Energy-efficient DC ECM Blower

- The leading energy efficient Class II Type B2 Biosafety Cabinet in the world with 70% energy savings compared to AC motor
- Stable airflow despite building voltage fluctuations and filter loading



## Advance Filtration System

- 10x Filtration efficiency of HEPA filter
- Creates ISO Class 3 work zone instead of Industry-standard ISO Class 5
- Same 10 years filter life and replacement cost as HEPA filters



## Dimmable LED

- Save energy and optimize work comfort




## Tray Support Beams

- Support work tray evenly for less vibration
- Cleaning holder to easily wipe the drain pan

## ISOCIDE™ Powder Coat

- Silver-ion impregnated powder coat
- Inhibits microbial growth to improve safety
- Prevents the plenum from becoming biohazard landfill



## Certification

	Performance	Air Quality	Filtration	Electrical Safety
Standards Compliance	NSF / ANSI 49, USA	ISO 14644.1, Class 3, Worldwide US Fed Std 209E, Class 1 USA JIS B9920, Class 3, Japan	EN-1822 (H14), Europe IEST-RP-CC001, USA	UL 61010-1 3rd Ed, USA CSA22.2, No.1010-192, Canada

## Dynamic Chamber™

- Blower plenum and side walls are surrounded by negative pressure
- Prevent contaminants from escaping outside


- Positive Pressure
- Negative Pressure




## Cabinet Filtration System

- Ambient air is pulled through the front grille to prevent contamination of the work surface and work product. The inflow does not mix with the clean air within the cabinet work zone.
- Ambient air is taken in through a pre-filter at the top of the cabinet, and passes through the downflow ULPA filter, entering the work zone as laminar flow. The uniform, non-turbulent air stream protects against cross contamination within and throughout the work area.
- Near the work surface, the downflow air stream splits with a portion moving toward the front air grille, and the remainder moving to the rear air grille. A small portion of the ULPA filtered downflow enters the intake perforations at the side capture zones at a higher velocity (small blue arrows).
- A combination of inflow and downflow air streams forms an air barrier that prevents contaminated room air from entering the work zone, and prevents work surface emissions from escaping the work zone. The downflow combined with the inflow air enters the common air plenum.
- All air in the common plenum is HEPA-filtered and exhausted via a dedicated ducting system to the external environment.
- Fail-safe system ensures that in case of exhaust failure, the cabinet's main fan automatically shuts down to ensure safety to the user

## The Performance Envelope Concept



- Nominal Airflow
- ▲ No Personnel / Product Protection
- Personnel / Product Protection
- Area of Personnel / Product Protection
- Area of no Personnel / Product Protection







Dynamic air barrier, where inflow and downflow converge  
Side capture zones

- ULPA-filtered air
- Unfiltered / potentially contaminated air
- Room air / Inflow air

## Engineering Drawing

314 mm / 360 mm / 579 mm / 797 mm  
(12.4" / 14.2" / 22.8" 31.4")



- VHP In Provision
- Pre-filter
- Exhaust Airflow Sensor
- USB Port
- Power Inlet
- Zero Voltage Relay Contact
- Centurion 7" Touch Screen Controllet

- Sash Glass
- Single-piece Stainless Steel Backwall
- Stainless Steel Work Tray
- Data Plate
- Dimmable LED Lamp
- IV Bar Retrofit Kit Provision

- Service Fixture Retrofit Kit provision
- Cable Port
- Stainless Steel Armrest
- Exhaust Collar
- Exhaust Filter
- DC ECM Blower
- Downflow Filter

- Downflow Sensor
- UV Lamp Provision
- Electrical outlet
- Drain Valve (optional)

Options and Accessories					
LB4-G4		3 feet	4 feet	5 feet	6 feet
Anti-blowback Valve 12 inches	EG Powder Coated	ABBV-12P 5170353			
	304 Stainless Steel	ABBV-12X 5170355			
UV Lamp		UV-15A 5170251	UV-30A 5170255		
IV Bar		IV-955 5170276	IV-1260 5170277	IV-1565 5170278	IV-1870 5170279
Electrical Outlet	Direct Mounted	EO-HC 5170035			
	GFCI	EO-GFCI 5170071			
Service Fixtures	EU SF-Gas-40 mm and Solenoid Valve	SF-1G40 5170002			
	EU SF-Vacuum-40 mm	SF-1V40 5170003			
	EU SF-Air-40 mm	SF-1A40 5170006			
	EU SF-Nitrogen-40 mm	SF-1N40 5170011			
	EU SF-Water-40 mm	SF-1W40 5170017			
	US SF-Universal-0 mm	SF-2U40 5170018			
	Copper Piping for SF	CU-Pipe 5170026			
Support Stand		STA-3A0 5131340	STA-4A0 5131341	STA-5A0 5131427	STA-6A0 5131389
Pipette Storage Shelf		STAINLESS STEEL PIPETTE STORAGE SHELF 5260327			
Arm Rest Padding		MEWREST 5170127			
Foot Rest		FT-REST 5170492			
Laboratory Chair		ME-LD-AR360 1150006			
IQOQ Protocol		IQOQ 9010179			



ABBV\_



UV\_A-L



IV\_



EO-H\_



EO-GFCI



SF-1\_



SF-2U\_



STA\_



Pipette Storage Shelf



MEWREST



FT-REST



ME-LD-AR360



IQOQ

## Class II Type B2 Biological Safety Cabinets (203 mm / 8" Opening)

TECHNICAL SPECIFICATIONS						
Labculture® Class II Type B2	Stainless Steel Side Walls	220-240 VAC, 50/60 Hz	LB2-3B8 G4 2011682	LB2-4B8 G4 2011684	LB2-5B8 G4 2011686	LB2-6B8 G4 2011688
		110-130 VAC, 50/60 Hz	LB2-3B9 G4 2011683	LB2-4B9 G4 2011685	LB2-5B9 G4 2011687	LB2-6B9 G4 2011689
Nominal Size			0.9 meter (3')	1.2 meter (4')	1.5 meter (5')	1.8 meter (6')
External Dimensions* (W x D x H)	Without Optional Base Stand		1115 x 867 x 1610 mm (43.9" x 34.1" x 63.4")	1420 x 867 x 1610 mm (55.9" x 34.1" x 63.4")	1725 x 867 x 1610 mm (67.9" x 34.1" x 63.4")	2030 x 867 x 1610 mm (79.9" x 34.1" x 63.4")
Internal Dimensions (W x D x H)			958 x 623 x 720 mm (37.7" x 24.5" x 28.3")	1263 x 623 x 720 mm (49.7" x 24.5" x 28.3")	1568 x 623 x 720 mm (61.7" x 24.5" x 28.3")	1873 x 623 x 720 mm (73.7" x 24.5" x 28.3")
Usable Work Area			0.45 m <sup>2</sup> (4.8 sq.ft.)	0.62 m <sup>2</sup> (6.7 sq.ft.)	0.76 m <sup>2</sup> (8.2 sq.ft.)	0.93 m <sup>2</sup> (10.0 sq.ft.)
Tested opening		203 mm (8")				
Maximum Sash Opening		534 mm (21")				
Average Airflow Velocity	Inflow	0.53 m/s (105 fpm)				
	Downflow	0.31 m/s (60 fpm)				
Airflow Volume	Inflow		376 m <sup>3</sup> /h (223 cfm)	493 m <sup>3</sup> /h (292 cfm)	608 m <sup>3</sup> /h (361 cfm)	724 m <sup>3</sup> /h (429 cfm)
	Downflow		628 m <sup>3</sup> /h (363 cfm)	822 m <sup>3</sup> /h (476 cfm)	1016 m <sup>3</sup> /h (588 cfm)	1210 m <sup>3</sup> /h (700 cfm)
	CBV Exhaust Air Volume**		1127 m <sup>3</sup> /h (658 cfm)	1476 m <sup>3</sup> /h (862 cfm)	1824 m <sup>3</sup> /h (1065 cfm)	2173 m <sup>3</sup> /h (1269 cfm)
	Minimum Exhaust Static Pressure		400 Pa / 1.6 in H <sub>2</sub> O	375 Pa / 1.5 in H <sub>2</sub> O	375 Pa / 1.5 in H <sub>2</sub> O	400 Pa / 1.6 in H <sub>2</sub> O
	CBV Exhaust Static Pressure**		575 Pa / 2.3 in H <sub>2</sub> O	550 Pa / 2.2 in H <sub>2</sub> O	550 Pa / 2.2 in H <sub>2</sub> O	575 Pa / 2.3 in H <sub>2</sub> O
Supply ULPA Filter Typical Efficiency		99.999% efficiency at 0.3-0.1 microns				
Exhaust HEPA Filter Typical Efficiency		≥99.99% at 0.3 microns				
Sound Emission (dBA)*	NSF / ANSI 49		57 dBA	58 dBA	59 dBA	60 dBA
	EN 12469		54 dBA	55 dBA	56 dBA	57 dBA
Light Intensity	LED Lamp Intensity	>1200 Lux (>93 foot-candles)				
	Optional UV Lamp	253.7 nm				
Electrical Requirements	Nominal Power Consumption		166 W	189 W	229 W	252 W
	Heat Load		566 BTU/Hr	645 BTU/Hr	781 BTU/Hr	860 BTU/Hr
	Full Load Amps	10 Ampere				
Cabinet Construction	Main Body	Electro-galvanized steel with white oven-baked epoxy-polyester Isocide™ antimicrobial powder-coated finish, 1.5 mm (0.06") / 16 gauge thick				
	Work Zone	Stainless steel Type 304 with No.4 finish, 1.5 mm (0.06") / 16 gauge thick				
Net Weight			279 Kg (615 lbs)	317 Kg (699 lbs)	359 Kg (791 lbs)	438 Kg (966 lbs)
Shipping Weight			318 Kg (703 lbs)	370 Kg (814 lbs)	402 Kg (886 lbs)	491 kg (1083 lbs)
Shipping Dimensions, Maximum (W x D x H)			1210 x 950 x 1950 mm (47.6" x 37.4" x 76.8")	1520 x 950 x 1950 mm (59.8" x 37.4" x 76.8")	1900 x 950 x 1950 mm (74.8" x 37.4" x 76.8")	2150 x 950 x 1950 mm (84.7" x 37.4" x 76.8")
Shipping Volume Dimensions (W x D x H)			2.24 m <sup>3</sup> (79.1 cu.ft.)	2.82 m <sup>3</sup> (99.6 cu.ft.)	3.52 m <sup>3</sup> (124.3 cu.ft.)	3.98 m <sup>3</sup> (140.6 cu.ft.)

Disclaimer: Technical Specifications may be subjected to further changes without further notice.

\*Electrical power consumption is an measurement of new unit with clean filter operated within nominal setpoint. Result per unit may vary.

\*\*This Concurrent Balance Value (CBV) Exhaust Volume (per Pitot Duct Traverse) and Static Pressure at cabinet exhaust connection should be used when sizing the HVAC exhaust and supply.