

Tube Chart for Rotor NU-RX500



MULTI-APPLICATION CENTRIFUGE



Centrifuge Model(s)	NU-C200V / NU-C300V	NU-C200R / NU-C300R / NU-C300RF
Max Speed RPM	4500	5100
Max Speed XG	4188	5350



Buckets NU-B500 ● ●




Lids NU-L500



Rotor Stand NU-RS2 ● ●

● **Cell Culture**
Accessory Package - C200G216

● **Blood Processing**
Accessory Package - C200G217

Tube Image													
Volume	3/5/7 ml	10 ml	15 ml	15 ml	50 ml	50 ml	100 ml	175 ml	225 ml	400 ml	625 ml		
Tube Shape	Round	Round	Round	Conical	Conical	Conical Skirted	Round	Conical	Conical	Round	Flat		
Max Diameter (mm)	13.5	17	17	17.5	30	30	47	62	62	80	93.4		
Max Height (mm)	117	117	117	124	124	124	133	130	130	125	133.5		
Insert Image													
Insert Part No.	NU-2T5 ●	NU-2T10 ●		NU-2T15C ●	NU-2T50CL	NU-2T50C ●	NU-2T50SL	NU-2T50S	NU-2T100	NU-2T200	NU-2T200	NU-2T400	NU-2T625
Qty Tubes/Insert	23	19	7	11	4	5	4	5	1	1	1	1	1
Qty Tubes/Rotor	92	76	28	44	16	20	16	20	4	4	4	4	4
Bio-containment lid NU-L500 Compatibility	Yes	Yes	Yes	Yes	Yes	No	Yes	No	Yes	Yes	No	Yes	No



Tube Chart for Rotor NU-RX250



MULTI-APPLICATION CENTRIFUGE



Rotor ● ●
NU-RX250



Buckets NU-B250 ● ●



Lids NU-L250



Rotor Stand NU-RS2 ● ●

- **Cell Culture Accessory Package** - C200G216
- **Blood Processing Accessory Package** - C200G217

Tube Image											
Volume	3/5/7 ml	10 ml	15 ml	15 ml	50 ml		50 ml		100 ml	175 ml	225 ml
Tube Shape	Round	Round	Round	Conical	Conical		Conical Skirted		Round	Conical	Conical
Max Diameter (mm)	13.5	17	17	17.5	30		30		47	62	62
Max Height (mm)	125	125	125	128	125		125		133	130	130
Insert Image											
Insert Part No.	NU-1T5 ●	NU-1T10 ●		NU-1T15C ●	NU-1T50CL	NU-2T50C ●	NU-1T50SL	NU-1T50S	NU-1T100	NU-1T200	NU-1T200
Qty Tubes/Insert	15	12	3	7	3	4	3	4	1	1	1
Qty Tubes/Rotor	60	48	12	28	12	16	12	16	4	4	4
Bio-containment lid NU-L500 Compatibility	Yes	Yes	Yes	Yes	Yes	No	Yes	No	Yes	Yes	No



Microplate and Angle Rotors



Rotor
NU-RX500



Buckets
NU-BXDW

Type of Tube	Part #	Description	Qty	Model: NU-C200V	Model: NU-C200R
96 Well Plate	NU-BXDW	Microplate Carrier (Set of 4)	16	4500 rpm 4188 xg	5100 rpm 5350 xg
Cell Culture Plate			16		
Deep Well			4		



NU-RA24-2



NU-RA48-2



NU-RA8-50



NU-RA30-15



NU-RA30-2



NU-RA16-5



Rotor NU-RHDW
Buckets NU-BH-DW
Bio-Lid NU-LH

Type of Tube	Part #	Description	Qty	Model: NU-C200V	Model: NU-C200R
Microtubes (1.5/2ml)	NU-RA24-2	Angular Rotor 24 x 1.5/2ml 45°	24	15300 RPM / 22380 xg	18000 RPM / 30980 xg
	NU-RA30-2	Angular Rotor 30 x 1.5/2ml 45°	30	15300 RPM / 25350 xg	16000 RPM / 27720 xg
	NU-RA48-2	Angular Rotor 48 x 1.5/2ml 45°	48	14000 RPM / 21720 xg	15500 RPM / 26620 xg
Microtubes (5ml)	NU-RA16-5	Angular Rotor 16 x 5ml 45°	16	13200 RPM / 21035 xg	13200 RPM / 21000 xg
Conical Tubes (15ml)	NU-RA30-15	Angular Rotor 30 x 15 37	30	4500 RPM / 3169 xg	4800 RPM / 3606 xg
Conical Tubes (50ml)	NU-RA8-50	Angular Rotor 8 x 50ml 25°	8	10700 RPM / 13000 xg	12100 RPM / 16660 xg

Type of Tube	Part #	Description	Qty	Model: NU-C200V	Model: NU-C200R
96 Well Plate	NU-RHDW	Microplate Rotor	8	4500 RPM / 3531 xg	5100 RPM / 4537 xg
Cell Culture Plate	NU-BH-DW	Microplate Buckets	8		
Deep Well	NU-LH	Microplate Lids	2		

