



# Product specification

**Aquatec** · Transfer system for microtomes of CUT series

Item No.: 10090056



## Specifications

Heated waterbath	up to +50 °C
Volume of heated waterbath	750 ml
Volume of pump reservoir	570 ml
Required liquid	demineralized or distilled water
Flow rate	0 - 500 ml / min.
Operating temperature range	+10 to +35 °C
Operating humidity	max. rel. 80 % non-condensing
Storage temperature range	+5 to +55 °C
Storage humidity	max. rel. 80 % non-condensing
Pollution degree	2
Overvoltage category	II
Power requirements	100 - 240 V / 50/60 Hz, 0.7 A + / -10 %
Weight	
basic device w/o transfer bridge, blade holder, incl. lid	5.0 kg
disposable blade holder AQ	2.4 kg
Dimensions (W x D x H) w/o transfer bridge, blade holder	
basic device mounted on microtome	300 x 235 x 128 mm
basic device unmounted on microtome	300 x 280 x 128 mm

## Standard components

Basic device incl. removable tray

---

Disposable blade holder AQ

---

Lightning

---

Waste tray

---

Holder for temperature sensor

---

Temperature sensor

---

2 x transfer inserts (1 x installed)

---

2 x suction filter (1 x installed)

---

Syringe

---

Drain screen

---

Manual

---

T: (763) 712-8717  
M: sales@gmi-inc.com