

VERDOT
SINGLE-USE
SYSTEMS

Worry Less.
Filter More.

Choose
FlexiPro™ TFF




VERDOT
SINGLE-USE
SYSTEMS

Choose Multi-purpose

GMP TFF processing from
clarification to formulation

**One System,
All You Need.**




VERDOT
SINGLE-USE
SYSTEMS

Choose Any Format

Compatible with **hollow fiber**
and **flat-sheet cassette**
filter formats

You Name It.
Flexipro Can Do It.



VERDOT
SINGLE-USE
SYSTEMS

Choose Configurability

Configurable flow kits with
advanced analytics for
better process control

**Your Process,
Your Way.**




VERDOT
SINGLE-USE
SYSTEMS

Choose Precision

Novel PCV combines precise backpressure and rapid response for **smooth TMP regulation**

**Less Fouling,
More Filtering.**



— Competitor — FlexiPro™ TFF

VERDOT
SINGLE-USE
SYSTEMS

Choose Intuitive

Directional grooves
streamline **mistake-proof**
Flow Kit installation

Ease of Use.
Ease of Mind.




VERDOT
SINGLE-USE
SYSTEMS

Choose Flexibility

Widest processing range
in a single **GMP** system
(1 to 3500 L/hr)

**Change the Flow Kit,
Not the System.**



Choose
FlexiPro™ TFF!

VERDOT

verdot-biotechnologies.com